

产品名称: Rab7 Rabbit Polyclonal Antibody
产品货号: APRab00106



产品概述 (Summary)

产品名称 (Production Name)	Rab7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB7A
别名 (Alternative Names)	RAB7; PRO2706
基因 ID (Gene ID)	7879
蛋白 ID (SwissProt ID)	P51149.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 23 kDa

研究背景 (Background)

Plays a role in the fusion of phagosomes with lysosomes. Plays important roles in microbial pathogen infection and

产品名称: Rab7 Rabbit Polyclonal Antibody
产品货号: AP Rab00106




survival, as well as in participating in the life cycle of viruses. Microbial pathogens possess survival strategies governed by RAB7A, sometimes by employing RAB7A function (e.g. Salmonella) and sometimes by excluding RAB7A function.


研究领域 (Research Area)

Tags & Cell Markers


图片 (Image Data)



Western blot analysis of RAB7 in (1) A375 lysates; (2) 3T3 lysates using Rab7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using RAB7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of Rab7 in HeLa using RAB7 antibody.

产品名称: Rab7 Rabbit Polyclonal Antibody
产品货号: APRab00106



注意事项 (Note)

For research use only .