

产品名称: MEK3/MEK6 Rabbit Polyclonal Antibody
产品货号: APRab00113



产品概述 (Summary)

产品名称 (Production Name)	MEK3/MEK6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP2K3/MAP2K6
别名 (Alternative Names)	MEK6; MKK6; MAPKK6; PRKMK6; SAPKK3; MAP2K6; MEK3; MAP kinase kinase 3; MAPKK3; MAPK/ERK kinase 3
基因 ID (Gene ID)	5606/5608
蛋白 ID (SwissProt ID)	P46734/P52564.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 39,37 kDa; Observed MW: 39,37 kDa

研究背景 (Background)

产品名称: MEK3/MEK6 Rabbit Polyclonal Antibody
产品货号: APRab00113




Dual specificity kinase. Is activated by cytokines and environmental stress in vivo. Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in the MAP kinase p38.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of MEK3/MEK6 in HeLa lysates using MEK3/MEK6 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using MEK3/MEK6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of MEK3/MEK6 in HeLa using MEK3/MEK6 antibody.

产品名称: MEK3/MEK6 Rabbit Polyclonal Antibody
产品货号: APRab00113



注意事项 (Note)

For research use only .