

产品名称: Trk Rabbit Polyclonal Antibody
产品货号: APRab00125



产品概述 (Summary)

产品名称 (Production Name)	Trk Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	NTRK1/NTRK2/NTRK3
别名 (Alternative Names)	gp140trk; GP145-TrkB; GP145-TrkC; MTC; NTRK1; NTRK2; NTRK3; p140-TrkA; TRKA; TRKB; TrkB tyrosine kinase; TRKC
基因 ID (Gene ID)	4914/4916/4915
蛋白 ID (SwissProt ID)	P04629/Q16288/Q16620.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 87 kDa; Observed MW: 140 kDa

研究背景 (Background)

产品名称: Trk Rabbit Polyclonal Antibody
产品货号: AP Rab00125




The family of Trk receptor tyrosine kinases consists of TrkA, TrkB, and TrkC. While the sequence of these family members is highly conserved, they are activated by different neurotrophins: TrkA by NGF, TrkB by BDNF or NT4, and TrkC by NT3. Neurotrophin signaling through these receptors regulates a number of physiological processes, such as cell survival, proliferation, neural development, and axon and dendrite growth and patterning.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)




Western blot analysis of TrkA+B+C in Human fetal brain lysates using Trk antibody.



Immunohistochemistry analysis of paraffin-embedded mouse brain using TrkA+B+C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Trk Rabbit Polyclonal Antibody
产品货号: APRab00125



Immunofluorescence analysis of Trk in Neuro2a using TrkA+B+C antibody.

注意事项 (Note)

For research use only .