

产品名称: JNK1/3 Rabbit Polyclonal Antibody
产品货号: APRab00131



产品概述 (Summary)

产品名称 (Production Name)	JNK1/3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAPK8/MAPK10
别名 (Alternative Names)	c-Jun N-terminal kinase 3; FLJ12099; FLJ33785; JNK1; JNK3; JNK3 alpha protein kinase; JNK3A; MAP kinase 10; MAP kinase p49 3F12; MAPK 10; MAPK10
基因 ID (Gene ID)	5599/5602
蛋白 ID (SwissProt ID)	P45983/P53779.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 48 kDa; Observed MW: 46,54 kDa

研究背景 (Background)

产品名称: JNK1/3 Rabbit Polyclonal Antibody
产品货号: APRab00131




The stress-activated protein kinase/Jun-amino-terminal kinase SAPK/JNK is potently and preferentially activated by a variety of environmental stresses including UV and gamma radiation, ceramides, inflammatory cytokines, and in some instances, by growth factors and GPCR agonists. As with the other MAPKs, the core signaling unit is composed of a MAPKKK, typically MEKK1-MEKK4, or by one of the mixed lineage kinases (MLKs), which phosphorylate and activate MKK4/7.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of JNK1/JNK3 in HeLa lysates using JNK1/3 antibody.

注意事项 (Note)

For research use only .