

产品名称: Peroxiredoxin 2 Rabbit Polyclonal Antibody
产品货号: AP Rab00135



产品概述 (Summary)

产品名称 (Production Name)	Peroxiredoxin 2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRDX2 PRDX2; NKEFB; TDPX1; Peroxiredoxin-2; Natural killer cell-enhancing factor B;
别名 (Alternative Names)	NKEF-B; PRP; Thiol-specific antioxidant protein; TSA; Thioredoxin peroxidase 1; Thioredoxin-dependent peroxide reductase 1
基因 ID (Gene ID)	7001
蛋白 ID (SwissProt ID)	P32119.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 22 kDa

产品名称: Peroxiredoxin 2 Rabbit Polyclonal Antibody
产品货号: AP Rab00135




研究背景 (Background)

This gene encodes a member of the peroxiredoxin family of antioxidant enzymes, which reduce hydrogen peroxide and alkyl hydroperoxides. The encoded protein may play an antioxidant protective role in cells, and may contribute to the antiviral activity of CD8(+) T-cells. This protein may have a proliferative effect and play a role in cancer development or progression.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Peroxiredoxin 2 in HEK293 lysates using Peroxiredoxin 2 antibody.

注意事项 (Note)

For research use only .