

产品名称: Fas Rabbit Polyclonal Antibody
产品货号: APRab00142



产品概述 (Summary)

产品名称 (Production Name)	Fas Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAS
别名 (Alternative Names)	FAS; ALPS1A; APO-1; APT1; CD95; FAS1; FASTM; TNFRSF6
基因 ID (Gene ID)	355
蛋白 ID (SwissProt ID)	P25445.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 48 kDa

研究背景 (Background)

Receptor for TNFSF6/FASLG. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-

产品名称: Fas Rabbit Polyclonal Antibody
产品货号: APRab00142




inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Fas in Ramos lysates using Fas antibody.

注意事项 (Note)

For research use only .