

产品名称: BMP11 Rabbit Polyclonal Antibody
产品货号: APRab00148



产品概述 (Summary)

产品名称 (Production Name)	BMP11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	GDF11
别名 (Alternative Names)	GDF11; BMP11; Bone morphogenetic protein 11; GDF-11; BMP-11
基因 ID (Gene ID)	10220
蛋白 ID (SwissProt ID)	O95390.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

Secreted signal that acts globally to specify positional identity along the anterior/posterior axis during development. Play

产品名称: BMP11 Rabbit Polyclonal Antibody
产品货号: APRab00148




critical roles in patterning both mesodermal and neural tissues and in establishing the skeletal pattern.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of BMP11 in (1) HeLa lysates; (2) C6 lysates using BMP11 antibody.

注意事项 (Note)

For research use only .