

产品名称: Sarcomeric alpha Actinin Rabbit Polyclonal Antibody
产品货号: APRab00151



产品概述 (Summary)

产品名称 (Production Name)	Sarcomeric alpha Actinin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTN2
别名 (Alternative Names)	Actin binding protein; Actinin alpha 2; ACTN2; Alpha actinin 2; CMD1AA
基因 ID (Gene ID)	88
蛋白 ID (SwissProt ID)	P35609.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 104 kDa; Observed MW: 104 kDa

研究背景 (Background)

产品名称: Sarcomeric alpha Actinin Rabbit Polyclonal Antibody
产品货号: AP Rab00151




F-actin cross-linking protein which is thought to anchor actin to a variety of intracellular structures. This is a bundling protein.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of ACTN2 in HepG2 lysates using Sarcomeric alpha Actinin antibody.

注意事项 (Note)

For research use only .