

产品名称: Ctip2 Rabbit Polyclonal Antibody  
产品货号: APRab00153



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | Ctip2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB                               |
| 种属反应性 (Reactivity)     | Human                            |

## 产品性能 (Performance)

|                     |                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                       |
| 修饰 (Modification)   | Unmodified                                                                                         |
| 同种型 (Isotype)       | IgG                                                                                                |
| 克隆 (Clonality)      | Polyclonal                                                                                         |
| 形式 (Form)           | Liquid                                                                                             |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Chromatography                                                                            |

## 免疫原信息 (Immunogen)

|                        |                                                                                                            |
|------------------------|------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | BCL11B                                                                                                     |
| 别名 (Alternative Names) | ATL1 alpha; ATL1 beta; ATL1 gamma; ATL1 delta; Bcl11b; BCL11B/TRDC fusion; CTIP2; hRIT1 alpha; hRit1; Rit1 |
| 基因 ID (Gene ID)        | 64919                                                                                                      |
| 蛋白 ID (SwissProt ID)   | Q9C0K0.                                                                                                    |

## 产品应用 (Application)

|                          |                                             |
|--------------------------|---------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000                             |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 96 kDa; Observed MW: 140 kDa |

## 研究背景 (Background)

**产品名称: Ctip2 Rabbit Polyclonal Antibody**  
**产品货号: APRab00153**




Tumor-suppressor protein involved in T-cell lymphomas. May function on the P53-signaling pathway. May be a key regulator of both differentiation and survival during thymocyte development. Repress transcription through direct, TFCOUP2-independent binding to a GC-rich response element.

## 研究领域 (Research Area)

Epigenetics and Nuclear Signaling

## 图片 (Image Data)



Western blot analysis of Ctip2 in Jurkat lysates using Ctip2 antibody.

## 注意事项 (Note)

For research use only .