

产品名称: Aryl Hydrocarbon Receptor Rabbit Polyclonal Antibody
产品货号: APRab00160



产品概述 (Summary)

产品名称 (Production Name)	Aryl Hydrocarbon Receptor Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	AHR
别名 (Alternative Names)	Ah receptor; AhR; Class E basic helix-loop-helix protein 76; bHLHe76; AHR
基因 ID (Gene ID)	196
蛋白 ID (SwissProt ID)	P35869.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 96 kDa; Observed MW: 105 kDa

研究背景 (Background)

**产品名称: Aryl Hydrocarbon Receptor Rabbit Polyclonal
Antibody**
产品货号: APRab00160




The aryl hydrocarbon receptor (AhR) is a ligand activated transcription factor involved in xenobiotic metabolism, cell cycle regulation, and development in response to both endogenous and environmental signals. Involved in cell-cycle regulation. Likely to play an important role in the development and maturation of many tissues.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of AhR in HEK293 lysates using Aryl Hydrocarbon Receptor antibody.

注意事项 (Note)

For research use only .