

产品名称: SDHA Rabbit Polyclonal Antibody  
产品货号: APRab00168



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | SDHA Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,IHC,ICC/IF,FC,IP             |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                 |

## 产品性能 (Performance)

|                     |                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                       |
| 修饰 (Modification)   | Unmodified                                                                                         |
| 同种型 (Isotype)       | IgG                                                                                                |
| 克隆 (Clonality)      | Polyclonal                                                                                         |
| 形式 (Form)           | Liquid                                                                                             |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Chromatography                                                                            |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                    |
|------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | SDHA                                                                                                                               |
| 别名 (Alternative Names) | SDHA; SDH2; SDHF; Succinate dehydrogenase [ubiquinone] flavoprotein subunit; mitochondrial; Flavoprotein subunit of complex II; Fp |
| 基因 ID (Gene ID)        | 6389                                                                                                                               |
| 蛋白 ID (SwissProt ID)   | P31040.                                                                                                                            |

## 产品应用 (Application)

|                          |                                                                             |
|--------------------------|-----------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 73 kDa; Observed MW: 73 kDa                                  |

## 研究背景 (Background)

**产品名称: SDHA Rabbit Polyclonal Antibody**  
**产品货号: APRab00168**




Flavoprotein (FP) subunit of succinate dehydrogenase (SDH) that is involved in complex II of the mitochondrial electron transport chain and is responsible for transferring electrons from succinate to ubiquinone.

## 研究领域 (Research Area)

Tags & Cell Markers

## 图片 (Image Data)



Western blot analysis of SDHA in (1) HeLa lysates; (2) mouse brain lysates; (3) rat brain lysates; using SDHA antibody.

## 注意事项 (Note)

For research use only .