

产品名称: NUR77 Rabbit Polyclonal Antibody
产品货号: APRab00195



产品概述 (Summary)

产品名称 (Production Name)	NUR77 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	NR4A1 NR4A1; GFRP1; HMR; NAK1; Nuclear receptor subfamily 4 group A member 1;
别名 (Alternative Names)	Early response protein NAK1; Nuclear hormone receptor NUR/77; Nur77; Orphan nuclear receptor HMR; Orphan nuclear receptor TR3; ST-59; Testicular receptor 3
基因 ID (Gene ID)	3164
蛋白 ID (SwissProt ID)	P22736.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 64 kDa; Observed MW: 64 kDa

产品名称: NUR77 Rabbit Polyclonal Antibody
产品货号: APRab00195




研究背景 (Background)

Orphan nuclear receptor. May act concomitantly with NURR1 in regulating the expression of delayed-early genes during liver regeneration. Binds the NGFI-B response element (NBRE) 5'-AAAAGGTCA-3' (By similarity). May inhibit NF-kappa-B transactivation of IL2.

研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of NUR77 in HepG2 lysates using NUR77 antibody.

注意事项 (Note)

For research use only .