

产品名称: Glucose Transporter GLUT3 Rabbit Polyclonal Antibody
产品货号: APRab00196



产品概述 (Summary)

产品名称 (Production Name)	Glucose Transporter GLUT3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC2A3
别名 (Alternative Names)	SLC2A3; GLUT3; Solute carrier family 2; facilitated glucose transporter member 3; Glucose transporter type 3; brain; GLUT-3
基因 ID (Gene ID)	6515
蛋白 ID (SwissProt ID)	P11169.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 54 kDa; Observed MW: 54 kDa

产品名称: Glucose Transporter GLUT3 Rabbit Polyclonal Antibody
产品货号: APRab00196




研究背景 (Background)

Facilitative glucose transporter. Probably a neuronal glucose transporter.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GLUT3 in (1) HepG2 lysates; (2) NIH/3T3 lysates using Glucose Transporter GLUT3 antibody.

注意事项 (Note)

For research use only .