

产品名称: alpha B Crystallin Rabbit Polyclonal Antibody
产品货号: AP Rab00197



产品概述 (Summary)

产品名称 (Production Name)	alpha B Crystallin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CRYAB CRYAB; CRYA2; Alpha-crystallin B chain; Alpha(B)-crystallin; Heat shock protein beta-5; HspB5; Renal carcinoma antigen NY-REN-27; Rosenthal fiber component
别名 (Alternative Names)	
基因 ID (Gene ID)	1410
蛋白 ID (SwissProt ID)	P02511.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 20 kDa; Observed MW: 20 kDa

产品名称: alpha B Crystallin Rabbit Polyclonal Antibody
产品货号: APRab00197




研究背景 (Background)

May contribute to the transparency and refractive index of the lens.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of CRYAB in Human fetal heart lysates using alpha B Crystallin antibody.

注意事项 (Note)

For research use only .