

**产品名称:** FRS2 Rabbit Polyclonal Antibody

**产品货号:** APRab00210

### 产品概述 (Summary)

|                   |                                                                                                    |
|-------------------|----------------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                                         |
| 宿主(Host)          | Rabbit                                                                                             |
| 应用(Application)   | WB,FC,IP                                                                                           |
| 种属反应(Reactivity)  | Human                                                                                              |
| 偶联物(Conjugation)  | Unconjugated                                                                                       |
| 修饰(Modification)  | Unmodified                                                                                         |
| 同种型(Isotype)      | IgG                                                                                                |
| 克隆(Clonality)     | Polyclonal                                                                                         |
| 剂型(Form)          | Liquid                                                                                             |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液(Buffer)      | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化 (Purification) | Affinity Chromatography                                                                            |

### 产品应用(Application)

|                       |                                            |
|-----------------------|--------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:1000,FC 1:50-1:100,IP 1:20-1:50 |
| 分子量(Molecular Weight) | Calculated MW: 57 kDa; Observed MW: 57 kDa |

### 抗原信息(Antigen Information)

|                       |                                                                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名(Gene Name)        | FRS2                                                                                                                                                |
| 别名(Alternative Names) | Fibroblast growth factor receptor substrate 2; FGFR substrate 2; FGFR-signaling adaptor<br>SNT; Suc1-associated neurotrophic factor target 1; SNT-1 |
| 基因 ID(Gene ID)        | 10818                                                                                                                                               |
| SwissProt ID          | Q8WU20                                                                                                                                              |
| 免疫原(Immunogen)        | Recombinant protein of human FRS2                                                                                                                   |

### 研究背景 (Background)

Adapter protein that links FGR and NGF receptors to downstream signaling pathways. Involved in the activation of MAP kinases. Modulates signaling via SHC1 by competing for a common binding site on NTRK1.


## 研究领域 (Research Area)

Signal Transduction

## 注意事项 (Note)

For research use only.

## 图片 (Image Data)



Western blot analysis of FRS2 in HeLa lysates using FRS2 antibody.

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>