

产品名称: Staufen Rabbit Polyclonal Antibody
产品货号: APRab00213



产品概述 (Summary)

产品名称 (Production Name)	Staufen Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAU1
别名 (Alternative Names)	STAU; PPP1R150
基因 ID (Gene ID)	6780
蛋白 ID (SwissProt ID)	O95793.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 63 kDa; Observed MW: 55 kDa

研究背景 (Background)

Binds double-stranded RNA (regardless of the sequence) and tubulin. May play a role in specific positioning of mRNAs at

产品名称: Staufen Rabbit Polyclonal Antibody
产品货号: APRab00213




given sites in the cell by cross-linking cytoskeletal and RNA components, and in stimulating their translation at the site.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Staufen in K562 lysates using Staufen antibody.

注意事项 (Note)

For research use only .