

产品名称: BNP Rabbit Polyclonal Antibody
产品货号: APRab00216



产品概述 (Summary)

产品名称 (Production Name)	BNP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPPB
别名 (Alternative Names)	NPPB; Natriuretic peptides B; Gamma-brain natriuretic peptide
基因 ID (Gene ID)	4879
蛋白 ID (SwissProt ID)	P16860.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 30 kDa

研究背景 (Background)

Brain natriuretic peptide (BNP) circulates in blood as a peptide hormone with natriuretic, vasodilatory and renin inhibitory

产品名称: BNP Rabbit Polyclonal Antibody
产品货号: APRab00216




properties. BNP is secreted predominantly by the left ventricular myocytes in response to volume expansion and pressure overload.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of BNP in Human plasma lysates using BNP antibody.

注意事项 (Note)

For research use only .