

产品名称: Adenosine A1 Receptor Rabbit Polyclonal Antibody
产品货号: APRab00224



产品概述 (Summary)

产品名称 (Production Name)	Adenosine A1 Receptor Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADORA1
别名 (Alternative Names)	ADORA1; Adenosine receptor A1
基因 ID (Gene ID)	134
蛋白 ID (SwissProt ID)	P30542.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 37 kDa

研究背景 (Background)

产品名称: Adenosine A1 Receptor Rabbit Polyclonal Antibody
产品货号: APRab00224




Receptor for adenosine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ADORA1 in A549 lysates using Adenosine A1 Receptor antibody.

注意事项 (Note)

For research use only .