

产品名称: Ran Rabbit Polyclonal Antibody
产品货号: AP Rab00230



产品概述 (Summary)

产品名称 (Production Name)	Ran Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAN RAN; ARA24; OK/SW-cl.81; GTP-binding nuclear protein Ran; Androgen receptor-associated protein 24; GTPase Ran; Ras-like protein TC4; Ras-related nuclear protein
别名 (Alternative Names)	
基因 ID (Gene ID)	5901
蛋白 ID (SwissProt ID)	P62826.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

产品名称: Ran Rabbit Polyclonal Antibody
产品货号: APRab00230




研究背景 (Background)

GTP-binding protein involved in nucleocytoplasmic transport. Required for the import of protein into the nucleus and also for RNA export. Involved in chromatin condensation and control of cell cycle (By similarity).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Ran in (1)HeLa lysates; (2)RAW264.7 lysates using Ran antibody.

注意事项 (Note)

For research use only .