

产品名称: Cyclophilin F Rabbit Polyclonal Antibody
产品货号: AP Rab00232



产品概述 (Summary)

产品名称 (Production Name)	Cyclophilin F Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPIF
别名 (Alternative Names)	PPIF; CYP3; Peptidyl-prolyl cis-trans isomerase F; mitochondrial; PPlase F; Cyclophilin F; Rotamase F
基因 ID (Gene ID)	10105
蛋白 ID (SwissProt ID)	P30405.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 18 kDa

研究背景 (Background)

产品名称: Cyclophilin F Rabbit Polyclonal Antibody
产品货号: APRab00232




PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of Cyclophilin F in HeLa lysates using Cyclophilin F antibody.

注意事项 (Note)

For research use only .