

产品名称: Fragilis Rabbit Polyclonal Antibody
产品货号: APRab00248



产品概述 (Summary)

产品名称 (Production Name)	Fragilis Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IFITM3
别名 (Alternative Names)	IFITM3; Interferon-induced transmembrane protein 3; Dispanin subfamily A member 2b; DSPA2b; Interferon-inducible protein 1-8U
基因 ID (Gene ID)	10410
蛋白 ID (SwissProt ID)	Q01628.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 15 kDa; Observed MW: 15 kDa

研究背景 (Background)

产品名称: Fragilis Rabbit Polyclonal Antibody
产品货号: APRab00248




IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens, namely influenza A H1N1 virus, West Nile virus (WNV), and dengue virus (WNV), by inhibiting the early step(s) of replication.

研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of Fragilis in HeLa lysates using Fragilis antibody.

注意事项 (Note)

For research use only .