

产品名称: PKC eta Rabbit Polyclonal Antibody
产品货号: APRab00249



产品概述 (Summary)

产品名称 (Production Name)	PKC eta Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRKCH
别名 (Alternative Names)	nPKC eta; PKC h; PKC L; PKCh; PKCL; Prkch; PRKCL; Protein kinase C eta
基因 ID (Gene ID)	5583
蛋白 ID (SwissProt ID)	P24723.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 78 kDa; Observed MW: 40,78 kDa

研究背景 (Background)

PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC also serves as the receptor

产品名称: PKC eta Rabbit Polyclonal Antibody
产品货号: AP Rab00249




for phorbol esters, a class of tumor promoters.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of PKC eta in MCF-7 lysates using PKC eta antibody.

注意事项 (Note)

For research use only .