

产品名称: OTX1/2 Rabbit Polyclonal Antibody
产品货号: APRab00250



产品概述 (Summary)

产品名称 (Production Name)	OTX1/2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ChIP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	OTX1/OTX2 Homeobox protein OTX1; Homeobox protein OTX2; Orthodenticle 1;
别名 (Alternative Names)	Orthodenticle 2; Orthodenticle homeobox 1; Orthodenticle homeobox 2; Orthodenticle1; Orthodenticle2; otx1; otx2
基因 ID (Gene ID)	5013/5015
蛋白 ID (SwissProt ID)	P32242/P32243.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ChIP 1:20
蛋白分子量 (Molecular Weight)	Calculated MW: 37 kDa; Observed MW: 32 kDa

产品名称: OTX1/2 Rabbit Polyclonal Antibody
产品货号: APRab00250




研究背景 (Background)

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence (BTS): 5'-TCTAATCCC-3'.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Otx1/2 in Y79 lysates using OTX1/2 antibody.

注意事项 (Note)

For research use only .