

产品名称: USP11 Rabbit Polyclonal Antibody
产品货号: APRab00254



产品概述 (Summary)

产品名称 (Production Name)	USP11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP11
别名 (Alternative Names)	UHX1; USP11
基因 ID (Gene ID)	8237
蛋白 ID (SwissProt ID)	P51784.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 110 kDa; Observed MW: 110 kDa

研究背景 (Background)

Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Inhibits the degradation of

产品名称: USP11 Rabbit Polyclonal Antibody
产品货号: APRab00254




target proteins by the proteasome. Plays a role in the regulation of pathways leading to NF-kappa-B activation. Plays a role in the regulation of DNA repair after double-stranded DNA breaks.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of USP11 in (1) HEK293 lysates; (2) mouse testis lysates using USP11 antibody.

注意事项 (Note)

For research use only .