

产品名称: FADD Rabbit Polyclonal Antibody
产品货号: APRab00258



产品概述 (Summary)

产品名称 (Production Name)	FADD Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	FADD FADD; MORT1; GIG3; Protein FADD; FAS-associated death domain protein;
别名 (Alternative Names)	FAS-associating death domain-containing protein; Growth-inhibiting gene 3 protein; Mediator of receptor induced toxicity
基因 ID (Gene ID)	8772
蛋白 ID (SwissProt ID)	Q13158.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 23 kDa; Observed MW: 28 kDa

产品名称: FADD Rabbit Polyclonal Antibody
产品货号: APRab00258




研究背景 (Background)

Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of FADD in A431 lysates using FADD antibody.

注意事项 (Note)

For research use only .