

产品名称: c-Rel Rabbit Polyclonal Antibody  
产品货号: APRab00261



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | c-Rel Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,IP                            |
| 种属反应性 (Reactivity)     | Human                            |

## 产品性能 (Performance)

|                     |                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                       |
| 修饰 (Modification)   | Unmodified                                                                                         |
| 同种型 (Isotype)       | IgG                                                                                                |
| 克隆 (Clonality)      | Polyclonal                                                                                         |
| 形式 (Form)           | Liquid                                                                                             |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.           |
| 储存溶液 (Buffer)       | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Chromatography                                                                            |

## 免疫原信息 (Immunogen)

|                        |         |
|------------------------|---------|
| 基因名 (Gene Name)        | REL     |
| 别名 (Alternative Names) | REL     |
| 基因 ID (Gene ID)        | 5966    |
| 蛋白 ID (SwissProt ID)   | Q04864. |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IP 1:20-1:50               |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 69 kDa; Observed MW: 69 kDa |

## 研究背景 (Background)

c-Rel contains an amino-terminal DNA-binding domain referred to as the REL homology domain (REH) and carboxy-

**产品名称: c-Rel Rabbit Polyclonal Antibody**  
**产品货号: AP Rab00261**




terminal transactivation domains. The c-Rel protein is typically inhibited in unstimulated cells by I $\kappa$ B $\alpha$  and I $\kappa$ B $\beta$ . c-Rel expression is highest in hematopoietic cells with extensive research studies demonstrating its role in immune cell function and pathogenesis of disease.

## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of c-Rel in Daudi lysates using c-Rel antibody.

## 注意事项 (Note)

For research use only .