

产品名称: CD166 Rabbit Polyclonal Antibody
产品货号: APRab00271



产品概述 (Summary)

产品名称 (Production Name)	CD166 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALCAM
别名 (Alternative Names)	ALCAM; MEMD; CD166 antigen; Activated leukocyte cell adhesion molecule; CD166
基因 ID (Gene ID)	214
蛋白 ID (SwissProt ID)	Q13740.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 65 kDa; Observed MW: 100 kDa

研究背景 (Background)

产品名称: CD166 Rabbit Polyclonal Antibody
产品货号: APRab00271




Cell adhesion molecule that binds to CD6. Involved in neurite extension by neurons via heterophilic and homophilic interactions. May play a role in the binding of T- and B-cells to activated leukocytes, as well as in interactions between cells of the nervous system.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD166 in Daudi lysates using CD166 antibody.

注意事项 (Note)

For research use only .