

产品名称: PSD93 Rabbit Polyclonal Antibody
产品货号: AP Rab00287



产品概述 (Summary)

产品名称 (Production Name)	PSD93 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	DLG2
别名 (Alternative Names)	DLG2; Disks large homolog 2; Channel-associated protein of synapse-110; Chapsyn-110; Postsynaptic density protein PSD-93
基因 ID (Gene ID)	1740
蛋白 ID (SwissProt ID)	Q15700.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 98 kDa; Observed MW: 110 kDa

研究背景 (Background)

产品名称: PSD93 Rabbit Polyclonal Antibody
产品货号: APRab00287




Required for perception of chronic pain through NMDA receptor signaling. Regulates surface expression of NMDA receptors in dorsal horn neurons of the spinal cord. Interacts with the cytoplasmic tail of NMDA receptor subunits as well as inward rectifying potassium channels.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of PSD93 in SH-SY5Y lysates using PSD93 antibody.

注意事项 (Note)

For research use only .