

产品名称: CDT2 Rabbit Polyclonal Antibody
产品货号: APRab00293



产品概述 (Summary)

产品名称 (Production Name)	CDT2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	DTL
别名 (Alternative Names)	CDW1; DCAF2; Dtl; L2DTL; RAMP
基因 ID (Gene ID)	51514
蛋白 ID (SwissProt ID)	Q9NZJ0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 79 kDa; Observed MW: 91 kDa

研究背景 (Background)

Substrate-specific adapter of a DCX (DDB1-CUL4-X-box) E3 ubiquitin-protein ligase complex required for cell cycle control,

产品名称: CDT2 Rabbit Polyclonal Antibody
产品货号: APRab00293




DNA damage response and translesion DNA synthesis. The DCX(DTL) complex, also named CRL4(CDT2) complex, mediates the polyubiquitination and subsequent degradation of CDT1 and CDKN1A/p21(CIP1).

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of CDT2 in 293T lysates using CDT2 antibody.

注意事项 (Note)

For research use only .