

产品名称: Glutamate Receptor 2/3 Rabbit Polyclonal Antibody
产品货号: APRab00295



产品概述 (Summary)

产品名称 (Production Name)	Glutamate Receptor 2/3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIA2/GRIA3
别名 (Alternative Names)	GluR K2; GluR K3; GLURB; GLURC; GRIA2; GRIA3
基因 ID (Gene ID)	2891/2892
蛋白 ID (SwissProt ID)	P42262/P42263.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 99 kDa

研究背景 (Background)

产品名称: Glutamate Receptor 2/3 Rabbit Polyclonal Antibody
产品货号: APRab00295




Ion channels activated by glutamate are typically divided into two classes. Those sensitive to N-methyl-D-aspartate (NMDA) are designated NMDA receptors (NMDAR) while those activated by α -amino-3-hydroxy-5-methyl-4-isoxalone propionic acid (AMPA) are known as AMPA receptors (AMPA).

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of GluR2+GluR3 in (1) Human fetal brain lysates; (2) mouse brain lysates using Glutamate Receptor 2/3 antibody.

注意事项 (Note)

For research use only .