

产品名称: Ikaros Rabbit Polyclonal Antibody
产品货号: APRab00299



产品概述 (Summary)

产品名称 (Production Name)	Ikaros Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	IKZF1 CLL associated antigen KW 6; DNA-binding protein Ikaros; hlk 1; hlk-1; Hs.54452; IK1; Ikaros (zinc finger protein); IKAROS; IKAROS family zinc finger 1 (Ikaros); Ikaros family zinc finger protein 1; Ikzf1; IKZF1_HUMAN; LYF1; Lymphoid transcription factor LyF-1; PRO0758; Zinc finger protein subfamily 1A 1 (Ikaros); zinc finger protein subfamily 1A 1; Zinc finger protein; subfamily 1A; member 1; ZNFN1A1.
别名 (Alternative Names)	
基因 ID (Gene ID)	10320
蛋白 ID (SwissProt ID)	Q13422.

产品应用 (Application)

产品名称: Ikaros Rabbit Polyclonal Antibody
产品货号: APRab00299



稀释比 (Dilution Ratio)

WB 1:500-1:1000, IHC 1:50-1:100, FC 1:50-1:100

蛋白分子量 (Molecular Weight)

Calculated MW: 58 kDa; Observed MW: 50-70 kDa


研究背景 (Background)

Transcription regulator of hematopoietic cell differentiation (PubMed:17934067). Binds gamma-satellite DNA (PubMed:17135265, PubMed:19141594). Plays a role in the development of lymphocytes, B- and T-cells.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Ikaros in MOLT4 lysates using Ikaros antibody.

注意事项 (Note)

For research use only .