

产品名称: CD134/OX40L Receptor Rabbit Polyclonal Antibody
产品货号: APRab00302



产品概述 (Summary)

产品名称 (Production Name)	CD134/OX40L Receptor Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNFRSF4
别名 (Alternative Names)	TNFRSF4; TXGP1L; Tumor necrosis factor receptor superfamily member 4; ACT35 antigen; OX40L receptor; TAX transcriptionally-activated glycoprotein 1 receptor; CD antigen CD134
基因 ID (Gene ID)	7293
蛋白 ID (SwissProt ID)	P43489.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa; Observed MW: 45 kDa

产品名称: CD134/OX40L Receptor Rabbit Polyclonal Antibody
产品货号: APRab00302




研究背景 (Background)

Receptor for TNFSF4/OX40L/GP34.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD134 in Human fetal spleen lysates using CD134/OX40L Receptor antibody.

注意事项 (Note)

For research use only .