

产品名称: NCK1 Rabbit Polyclonal Antibody
产品货号: APRab00306



产品概述 (Summary)

产品名称 (Production Name)	NCK1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCK1
别名 (Alternative Names)	NCK; nck-1; NCKalpha
基因 ID (Gene ID)	4690
蛋白 ID (SwissProt ID)	P16333.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa; Observed MW: 47 kDa

研究背景 (Background)

Adapter protein which associates with tyrosine-phosphorylated growth factor receptors or their cellular substrates.

产品名称: NCK1 Rabbit Polyclonal Antibody
产品货号: APRab00306




Maintains low levels of EIF2S1 phosphorylation by promoting its dephosphorylation by PP1.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Nck in HeLa lysates using NCK1 antibody.

注意事项 (Note)

For research use only .