

产品名称: DNA Nucleotidylexotransferase Rabbit Polyclonal Antibody
产品货号: APRab00309



产品概述 (Summary)

产品名称 (Production Name)	DNA Nucleotidylexotransferase Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNTT
别名 (Alternative Names)	DNA nucleotidylexotransferase; DNTT; TDT; Terminal addition enzyme; Terminal deoxynucleotidyltransferase; Terminal transferase
基因 ID (Gene ID)	1791
蛋白 ID (SwissProt ID)	P04053.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 59 kDa; Observed MW: 59 kDa

产品名称: DNA Nucleotidylexotransferase Rabbit Polyclonal Antibody
产品货号: AP Rab00309




研究背景 (Background)

Template-independent DNA polymerase which catalyzes the random addition of deoxynucleoside 5'-triphosphate to the 3'-end of a DNA initiator. One of the in vivo functions of this enzyme is the addition of nucleotides at the junction (N region) of rearranged Ig heavy chain and T-cell receptor gene segments during the maturation of B- and T-cells.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of TDT in Jurkat lysates using DNA Nucleotidylexotransferase antibody.

注意事项 (Note)

For research use only .