

产品名称: Cytochrome P450 Reductase Rabbit Polyclonal Antibody
产品货号: APRab00320



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 Reductase Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	POR
别名 (Alternative Names)	CPR; CYPOR; P450R; por
基因 ID (Gene ID)	5447
蛋白 ID (SwissProt ID)	P16435.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 77 kDa; Observed MW: 77 kDa

研究背景 (Background)

产品名称: Cytochrome P450 Reductase Rabbit Polyclonal Antibody
产品货号: AP Rab00320




This enzyme is required for electron transfer from NADP to cytochrome P450 in microsomes. It can also provide electron transfer to heme oxygenase and cytochrome B5.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Cytochrome P450 Reductase in A431 lysates using Cytochrome P450 Reductase antibody.

注意事项 (Note)

For research use only .