

产品名称: Geminin Rabbit Polyclonal Antibody
产品货号: APRab00322



产品概述 (Summary)

产品名称 (Production Name)	Geminin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	GMNN
别名 (Alternative Names)	Gem; Geminin; GMNN
基因 ID (Gene ID)	51053
蛋白 ID (SwissProt ID)	O75496.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 30 kDa

研究背景 (Background)

Inhibits DNA replication by preventing the incorporation of MCM complex into prereplication complex (pre-RC). It is

产品名称: Geminin Rabbit Polyclonal Antibody
产品货号: APRab00322




degraded during the mitotic phase of the cell cycle. Its destruction at the metaphase-anaphase transition permits replication in the succeeding cell cycle.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Geminin in HeLa lysates using Geminin antibody.

注意事项 (Note)

For research use only .