

产品名称: LIS1 Rabbit Polyclonal Antibody
产品货号: APRab00327



产品概述 (Summary)

产品名称 (Production Name)	LIS1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	PAFAH1B1
别名 (Alternative Names)	MDS; LIS1; LIS2; MDCR; NudF; PAFAH
基因 ID (Gene ID)	5048
蛋白 ID (SwissProt ID)	P43034.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 47 kDa; Observed MW: 47 kDa

研究背景 (Background)

Required for proper activation of Rho GTPases and actin polymerization at the leading edge of locomoting cerebellar

产品名称: LIS1 Rabbit Polyclonal Antibody
产品货号: APRab00327




neurons and postmigratory hippocampal neurons in response to calcium influx triggered via NMDA receptors.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of LIS1 in SH-SY5Y lysates using LIS1 antibody.

注意事项 (Note)

For research use only .