

产品名称: SLC9A1 Rabbit Polyclonal Antibody
产品货号: APRab00334



产品概述 (Summary)

产品名称 (Production Name)	SLC9A1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC9A1
别名 (Alternative Names)	APNH; APNH1; NHE1; SLC9A1
基因 ID (Gene ID)	6548
蛋白 ID (SwissProt ID)	P19634.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa; Observed MW: 130 kDa

研究背景 (Background)

Involved in pH regulation to eliminate acids generated by active metabolism or to counter adverse environmental

产品名称: SLC9A1 Rabbit Polyclonal Antibody
产品货号: APRab00334




conditions. Major proton extruding system driven by the inward sodium ion chemical gradient. Plays an important role in signal transduction.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Sodium / Hydrogen Exchanger 1 in MOLT4 lysates using SLC9A1 antibody.

注意事项 (Note)

For research use only .