

产品名称: MCM3 Rabbit Polyclonal Antibody
产品货号: APRab00336



产品概述 (Summary)

产品名称 (Production Name)	MCM3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCM3
别名 (Alternative Names)	MCM3; HCC5; P1-MCM3; P1.h; RLFB
基因 ID (Gene ID)	4172
蛋白 ID (SwissProt ID)	P25205.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 91 kDa; Observed MW: 100 kDa

研究背景 (Background)

Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once

产品名称: MCM3 Rabbit Polyclonal Antibody
产品货号: APRab00336




per cell cycle' DNA replication initiation and elongation in eukaryotic cells.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of MCM3 in (1) K562 lysates; (2) RAW264.7 lysates using MCM3 antibody.

注意事项 (Note)

For research use only .