

产品名称: Coilin Rabbit Polyclonal Antibody
产品货号: APRab00340



产品概述 (Summary)

产品名称 (Production Name)	Coilin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	COIL
别名 (Alternative Names)	CLN80; COIL; Coilin; Coilin p80
基因 ID (Gene ID)	8161
蛋白 ID (SwissProt ID)	P38432.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 63 kDa; Observed MW: 80 kDa

研究背景 (Background)

Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of

产品名称: Coilin Rabbit Polyclonal Antibody
产品货号: APRab00340




nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of Coilin in Jurkat lysates using Coilin antibody.

注意事项 (Note)

For research use only .