

产品名称: E2F4 Rabbit Polyclonal Antibody
产品货号: APRab00342



产品概述 (Summary)

产品名称 (Production Name)	E2F4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	E2F4
别名 (Alternative Names)	E2F4; Transcription factor E2F4; E2F-4
基因 ID (Gene ID)	1874
蛋白 ID (SwissProt ID)	Q16254.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 44 kDa; Observed MW: 60 kDa

研究背景 (Background)

Transcription activator that binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3'

产品名称: E2F4 Rabbit Polyclonal Antibody
产品货号: APRab00342




found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of E2F4 in Raji lysates using E2F4 antibody.

注意事项 (Note)

For research use only .