

产品名称: STAT6 Rabbit Polyclonal Antibody
产品货号: APRab00349



产品概述 (Summary)

产品名称 (Production Name)	STAT6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAT6
别名 (Alternative Names)	STAT6; Signal transducer and activator of transcription 6; IL-4 Stat
基因 ID (Gene ID)	6778
蛋白 ID (SwissProt ID)	P42226.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 94 kDa; Observed MW: 110 kDa

研究背景 (Background)

STAT6 transcription factor of the STAT family. Plays a central role in IL4-mediated biological responses. Induces the

产品名称: STAT6 Rabbit Polyclonal Antibody
产品货号: APRab00349




expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of STAT6 in Raw264.7 lysates using STAT6 antibody.

注意事项 (Note)

For research use only .