

产品名称: CCR8 Rabbit Polyclonal Antibody
产品货号: APRab00351



产品概述 (Summary)

产品名称 (Production Name)	CCR8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC
种属反应性 (Reactivity)	Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCR8
别名 (Alternative Names)	CKRL1; CMKBR8; CMKBRL2
基因 ID (Gene ID)	1237
蛋白 ID (SwissProt ID)	P51685.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

Receptor for the chemokine CCL1/SCYA1/I-309. May regulate monocyte chemotaxis and thymic cell line apoptosis.

产品名称: CCR8 Rabbit Polyclonal Antibody
产品货号: AP Rab00351




Alternative coreceptor with CD4 for HIV-1 infection.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CCR8 in rat large intestine lysates using CCR8 antibody.

注意事项 (Note)

For research use only .