

产品名称: IKK epsilon Rabbit Polyclonal Antibody
产品货号: APRab00354



产品概述 (Summary)

产品名称 (Production Name)	IKK epsilon Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IKBKE IKBKE; IKKE; IKKI; KIAA0151; Inhibitor of nuclear factor kappa-B kinase subunit epsilon; I-kappa-B kinase epsilon; IKK-E; IKK-epsilon; IkBKE; Inducible I kappa-B kinase; IKK-i
别名 (Alternative Names)	
基因 ID (Gene ID)	9641
蛋白 ID (SwissProt ID)	Q14164.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 80 kDa; Observed MW: 80 kDa

产品名称: IKK epsilon Rabbit Polyclonal Antibody
产品货号: APRab00354




研究背景 (Background)

Phosphorylates inhibitors of NF-kappa-B thus leading to the dissociation of the inhibitor/NF-kappa-B complex and ultimately the degradation of the inhibitor. May play a special role in the immune response. Protects cells against DNA damage-induced cell death.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of IKKi/IKKe in Hela, Raji, Ramos, MCF-7 lysates using IKK epsilon antibody.

注意事项 (Note)

For research use only .