

产品名称: EMSY Rabbit Polyclonal Antibody
产品货号: APRab00355



产品概述 (Summary)

产品名称 (Production Name)	EMSY Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EMSY
别名 (Alternative Names)	BRCA2-interacting transcriptional repressor EMSY
基因 ID (Gene ID)	56946
蛋白 ID (SwissProt ID)	Q7Z589.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 141 kDa; Observed MW: 141 kDa

研究背景 (Background)

Regulator which is able to repress transcription, possibly via its interaction with a multiprotein chromatin remodeling

产品名称: EMSY Rabbit Polyclonal Antibody
产品货号: APRab00355




complex that modifies the chromatin.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of EMSY in HL-60, U2OS, C6 lysates using EMSY antibody.

注意事项 (Note)

For research use only .