

产品名称: Lipin 1 Rabbit Polyclonal Antibody
产品货号: APRab00356



产品概述 (Summary)

产品名称 (Production Name)	Lipin 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPIN1
别名 (Alternative Names)	Lipin-1; Lpin1; PAP1
基因 ID (Gene ID)	23175
蛋白 ID (SwissProt ID)	Q14693.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 130 kDa

研究背景 (Background)

Plays important roles in controlling the metabolism of fatty acids at different levels. Acts as a magnesium-dependent

产品名称: Lipin 1 Rabbit Polyclonal Antibody
产品货号: AP Rab00356




phosphatidate phosphatase enzyme which catalyzes the conversion of phosphatidic acid to diacylglycerol during triglyceride, phosphatidylcholine and phosphatidylethanolamine biosynthesis in the reticulum endoplasmic membrane.

研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Lipin 1 in Jurkat lysates using Lipin 1 antibody.

注意事项 (Note)

For research use only .