

产品名称: GIV Rabbit Polyclonal Antibody
产品货号: APRab00358



产品概述 (Summary)

产品名称 (Production Name)	GIV Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, FC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCDC88A
别名 (Alternative Names)	Akt phosphorylation enhancer; APE; Coiled-coil domain-containing protein 88A; G alpha-interacting vesicle-associated protein; GIV; HkRP1
基因 ID (Gene ID)	55704
蛋白 ID (SwissProt ID)	Q3V6T2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, FC 1:50-1:100, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 216 kDa; Observed MW: 216 kDa

研究背景 (Background)

产品名称: GIV Rabbit Polyclonal Antibody
产品货号: APRab00358




Bifunctional modulator of guanine nucleotide-binding proteins (G proteins).


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of GIV in mouse brain, Raw264.7 lysates using GIV antibody.



Western blot analysis of Girdin in Hela, A549, HL-60, U2OS, C6 lysates using Girdin antibody.

注意事项 (Note)

For research use only .