

产品名称: CD79B Rabbit Polyclonal Antibody
产品货号: APRab00359



产品概述 (Summary)

产品名称 (Production Name)	CD79B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD79B
别名 (Alternative Names)	CD79B; AGM6; B29; IGB
基因 ID (Gene ID)	974
蛋白 ID (SwissProt ID)	P40259.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 30-50 kDa

研究背景 (Background)

Required in cooperation with CD79A for initiation of the signal transduction cascade activated by the B-cell antigen

产品名称: CD79B Rabbit Polyclonal Antibody
产品货号: APRab00359




receptor complex (BCR) which leads to internalization of the complex, trafficking to late endosomes and antigen presentation.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD79B in Ramos, rat spleen lysates using CD79B antibody.

注意事项 (Note)

For research use only .