

产品名称: KA1 Rabbit Polyclonal Antibody
产品货号: APRab00360



产品概述 (Summary)

产品名称 (Production Name)	KA1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIK4
别名 (Alternative Names)	EAA1; GluK4; GRIK; Grik4; KA1
基因 ID (Gene ID)	2900
蛋白 ID (SwissProt ID)	Q16099.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 107 kDa; Observed MW: 107 kDa

研究背景 (Background)

Receptor for glutamate. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system.

产品名称: KA1 Rabbit Polyclonal Antibody
产品货号: APRab00360




The postsynaptic actions of Glu are mediated by a variety of receptors that are named according to their selective agonists.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of KA1 in mouse brain lysates using KA1 antibody.

注意事项 (Note)

For research use only .