

产品名称: LI Cadherin Rabbit Polyclonal Antibody
产品货号: APRab00363



产品概述 (Summary)

产品名称 (Production Name)	LI Cadherin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDH17
别名 (Alternative Names)	Cadherin 17; Cdh17; HPT 1; LI-cadherin; Liver Cadherin; liver inTine cadherin
基因 ID (Gene ID)	1015
蛋白 ID (SwissProt ID)	Q12864.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa; Observed MW: 120 kDa

研究背景 (Background)

Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic

产品名称: LI Cadherin Rabbit Polyclonal Antibody
产品货号: APRab00363




manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. LI-cadherin may have a role in the morphological organization of liver and intestine. Involved in intestinal peptide transport.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of LI Cadherin in HCT116, SW480 lysates using LI Cadherin antibody.

注意事项 (Note)

For research use only .